



WHEN TECHNOLOGY opens *the way* !



The new vascular
ACCESS EQUATION
4 x 9F = 20F

www.alseal.net

ALSeal proposes you a **global solution** of medical devices
CE marked for the management of endovascular percutaneous access.



HQS : Ease of use

- Single hand valve operation
- Secure locking
- Guidewire centering device included
- Side port included.

Available **from 14 to 22F** with a useful length of 140 and 190 mm.

XCath

Allows multiple access.

HQS Active Valve introducer

To manage the vascular access with an efficient sealing

Minimal blood loss for aortic endoluminal procedure.
Sealing on guidewires through endovascular tools up to 22F.



XCath

To convert the HQS Active Valve into multiple access

X-Cath device has been specifically designed to adapt HQS Introducer into multiple accesses.
X-Cath device ensures the sealing all over the procedure, using simultaneously up to 4 tools with the same HQS introducer.



Product sizes and guide


When technology opens the way !

HQS Active Valve Introducer

		HQS Active Valve Introducer (useful Ø)				
		14F	16F	18F	20F	22F
HQS Active Valve Introducer (length)	140 mm	✓	✓	✓	✓	✓
	190 mm			✓	✓	✓

HQS Active Valve Introducer product range reference: 2064-HQS

XCath 2064-XCATH-4 4 x 9F

- All of the products are supplied separately, sterile, for single use only.
- HQS Active Valve Introducer and XCath are unitary packaged
- Medical devices marked  2797.

Medical device Class IIa



Temis Santé - Bioparc

2 rue du Professeur Milleret
25000 Besançon - France

Phone: +33 (0)3 81 61 69 93

Fax: +33 (0)3 81 53 47 65

@ commercial@alcis.net
jp.barbault@alcis.net

www.alseal.net

